


Reduce your operation time

with the Tilly Medical

Bead Injector



Bead marking made quick and easy

With the new Bead Injector the operation time required for implanting tantalum beads is significantly reduced. Using the Bead Injector with its revolving magazine for tantalum beads, you will have the next bead in the needle within a second!

Easy Handling

You can use the loaded injector as a bradawl to position the bead. The Bead Injector is a strong, durable precision instrument with a needle support sleeve, but being made of titanium it is nevertheless lightweight.

The magazine is easy to load with prepacked 27-bead containers.

Technical data:

- Dimensions: 70 x 28 x 42 mm
- Weight: 90 gr
- Material: Titan
- Magazine capacity: 27 beads
- Ergonomic
- Easy to assemble/disassemble

The Bead Injector has been developed in co-operation with the Department of Orthopaedics, Lund University Hospital, Sweden.

For further information please contact:

Tilly Medical Products AB

Ideon Science & Technology Park, 223 70 Lund, Sweden

Telephone: +46 46 13 10 10, Fax: +46 46 13 10 33

Email: mail@tmed.se WWW: <http://www.tmed.se>